

# BODYFILL 60 iREJU

Cross-linked Hyaluronic Acid  
HA 24mg/mL



**The highly purified** viscoelastic body filler iREJU BODY HA Filler is the best solution for replenishing body volume and professional body contouring. iREJU BODY HA FILLER **restores volume to the body** in areas like breasts, lower legs, and buttocks. It is also efficient in combating cellulite and uneven skin surface due to liposuction procedures.

## 3C Technology of iREJU

### Conservative Cross-linking system

The 14-day process, low-temperature crosslinking process slowly mixes the crosslinked molecular weight and free HA to ensure uniform texture with 1-2 micron particles and preserves viscoelastic properties that allow for better tissue integration.

### Certified safe raw material

We use high-purity HA whose safety has been fully verified worldwide.

### Conservative Cross-linking system

Minimize differences between production batches through strict quality control

Hyaluronic acid in iREJU Body HA Filler has the highest level of purification for maximum safety, leaving minimal residue of cross-linking agent, endotoxins, and proteins.


## Benefits

- iREJU BODY HA Filler cross-linking that ensures a longer effect
- Contains highly-purified stabilized hyaluronic acid
- Has higher viscosity that prevents deformation
- Smooths the skin and combats cellulite
- Suitable for patients, who are not intended to go through surgery, but desire to enlarge volume in areas like breasts and buttocks
- Hyaluronic acid gel does not contain animal components, eliminating the risk of allergic reactions
- General anesthesia is not required
- No incisions or scars
- Quick recovery period (compared with plastic surgery)

## Indications

- Asymmetrical breasts and buttocks
- Uneven skin surface, including cellulite
- Volume loss in breasts and buttocks
- Contour correction of genitals (enlargement of labia majora, penis)
- Subcutaneous tissue loss due to liposuction
- Funnel chest
- Asymmetric calves
- Soft tissue defects related to surgery and trauma

## Unique Technology



- High elasticity
- High quality
- Safety
- Natural volume

## Characteristics

The experimental group showed better physical properties than the control group.

- High cross-linked HA gel content compared to similar products
- High viscoelasticity of cross-linked HA gel suitable for body filler

Products	HA Content (mg/ml)	G' (Pa) Storage modulus	$\eta^*$ (Pa.S) Complex viscosity	G'' (Pa) Loss modulus	Tan Delta (G'/G'')
iREJU	24	470.6	430.0	92.8	0.20
Control A	20	395.4	505.5	43.7	0.11
Control B	20	212.4	435.3	41.1	0.19
Control C	20	163.8	205.6	28.2	0.17



## NOTICE

The majority of our products are designed for professional use by doctors and specialists, by purchasing said product(s) you agree to supply on request formal proof of medical/professional training. Said product(s) are not intended for personal or home use. Please note this product has contraindications for use. It is recommended that all Healthcare professionals familiarise themselves with these contraindications and precautions by reading the instructions provided.



sbodyline.com



Restores Volume  
to the breasts  
without surgery



Volumizes the  
buttocks along  
with minimally  
invasive method



## Body Contour Filler "iREJU BODYFILL 60"



**DFKBIOLAB INC.**

**ADDRESS.** #522, 272 Sunhwagung-ro, Namyangju-si, Gyeonggi-do, 12106, Republic of Korea

**TEL.** +82-2-6225-7988 (7985,7984) **E-mail.** dfkbiolab@sbodyline.com **Homepage.** www.sbodyline.com